

	No	Type		No	Type		No	Type
AB Dental			Bio3			Cowellmedi		
Narrow Platform	55021	11	Advanced	55067	59	Inno Int Octa	55014	4
Standard Platform	55049	41	Progressive	55067	59	Inno Sub Hex	55037	29
ADIN			BioHorizons			C-Tech		
Touareg S/OS	55034	26	Laser-Lok 3.0 Grey	55034	26	EL	55067	59
Touareg Closefit UNP 2.75	55059	51	Internal Laser-Lok 3.5 Yellow	55049	41			
Touareg Closefit NP 3.0	55059	51	Internal Laser-Lok 4.5 Green	55051	43	Dentack		
Touareg Closefit RP 3.5	55046	38	Internal Laser-Lok 5.7 Blue	55051	43	Internal Hex	55056	48
Touareg CloseFit WP 4.3, 5.0	55067	59	External 4.0 Platform	55011	1			
Swell 3.3	55034	26				Dentalis		
			Biomet/3i			Regular	55051	43
AIRES			External Hex	55011	1			
3.3	55049	41	Certain internal Purple	55036	28	Dental Ratio		
			Certain internal Purple Platform 3.4	55036	28	Tissue Level RP, WP	55014	4
			Certain internal Blue, Yellow, Green	55066	58	Bone Level NC	55061	53
Alpha-Bio			T3 Purple	55036	28	Bone Level RC	55062	54
Internal Hex	55034	26	T3 Blue, Yellow, Green	55066	58			
Conical Hex 3.2	55065	57				Dental Tech		
Conical Hex 3.75-5.0	55051	43	Bionnovation			Implassic FT3 3.25	55039	31
			Biomorse	55023	14	Implassic FT3 3.75	55049	41
Alfa Gate			Biotech Dental			Dentegris		
Bioactive, All diam	55034	26	Contact, Contact-S 3.6 - 4.8	55045	37	3.3	55034	26
MAX, All diam	55034	26				3.75	55057	49
M+ NP 3.5	55044	36	BlueSkyBio			4.5	55049	41
M+ RP 4.3 - 5.0	55037	29	BIO Max NP Platform, All	55046	38	3.8/4.5	55049	41
Porous, All diam	55034	26	BIO Max RP Platform 4.3-5.0	55037	29	5.5	55020	10
			BIO Int Hex 3.5 Platform 3.7-4.1	55049	41			
Alpha Dent			BIO Int Hex 3.5 Platform 4.7	55051	43	Dentis		
Active Conus	55067	59	BIO Int Hex 4.5 Platform 4.7-7.0	55051	43	Cleanlant S-Clean 3.7 - 4.8	55020	10
Classic Conus	55067	59	BIO One-stage Narrow, Regular, Wide	55014	4	S-Clean OneQ Narrow 3.3	55064	56
Active Bio	55067	59	BIO One-stage Short Regular & Wide	55014	4			
Slim	55039	31	BIO Quattro Narrow 3.3	55061	53	Dentium		
			BIO Quattro Regular 4.1 - 5.6	55062	54	Superline	55037	29
Anker			BIO Conus 12, 3.0 Platform 3.0	55064	56	NR Line	55066	58
SB III 3.5	55023	14	BIO Conus 12, 3.5/4.0 Platform 3.5-4.0	55046	38	Simple Line II GH 1,6	55014	4
SB III 4.0-5.0	55020	10	BIO Conus 12, 4.5 x 6mm	55046	38			
			BIO Conus 12, 4.5/5.0 Platform 4.5 - 5.0	55017	7	DEVA		
			BIO Conus 12, 5.0 x 6mm	55046	38	Bone Level 3.5	55064	56
Ankylos			BIO Trilobe 3.5	55016	6	Bone Level 4.5	55051	43
C/X	55025	16	BIO Trilobe 4.3-6.0	55017	7			
			Bredent			DIO		
Anthogyr			SKY	55025	16	SM Submerged 3.8	55049	41
Axiom Bone Level 3.4 - 5.2	55065	57	copaSKY	55027	19	SM Submerged 5.0	55020	10
Axiom Tri-lobe	55065	57	mini2SKY	55033	25	UF II Narrow	55064	56
Axiom Tissue Level Narrow, Regular	55050	42	fast & fixed Abutment	55083	75	UF II Regular	55020	10
			BTI			Easy Implant		
Argon			Interna Universal	55051	43	Ocean	55046	38
K3pro Rapid 3.0 - 3.5	55026	17	Interna Wide	55051	43	Lilas	55037	29
K3pro Rapid Platf 3mm	55019	9	BTK			Edierre		
K3pro Rapid Platf 3mm Short	55044	36	External Hex EN	55038	30	3-lobe	55021	11
			External Hex ER	55011	1	4-lobe	55025	16
Astra Tech			External Hex EW	55011	1			
Osseospeed EV 3.0	55058	50	Internal Hex IR	55049	41	Euroteknika		
Osseospeed EV 3.6	55046	38	Internal Hex IM	55051	43	Internal Platform 3.0	55064	56
Osseospeed EV 4.2	55057	49	Internal Hex IW	55051	43	Internal Platform 3.5	55081	73
Osseospeed EV 4.8	55055	47	Hexagonal Morse Taper KR	55081	73	Naturactis/Naturall+ 3.0	55064	56
Osseospeed EV 5.4	55048	40	Hexagonal Morse Taper KW	55056	48	Naturactis 3.5-5.0	55081	73
Osseospeed TX 3.0 S	55064	56	Octagonal Morse Taper SR	55014	4	Uneva	55011	1
Osseospeed TX 3.5-4.0 S	55046	38	Camlog			galimplant		
Osseospeed TX 4.5-5.0 S	55017	7	Screw-Line, Root-Line 3.3	55036	28	iPx 4012	55046	38
			Screw-Line, Root-Line 3.8, 4.3	55031	23			
Avinent			Screw-Line, Root-Line 5.0, 6.0	55032	24	GC Tech		
External, Coral & Ocean 3.5	55038	30	iSy	55065	57	Aadva Narrow	55081	73
External, Coral & Ocean 4.1	55011	1	Ceralog			Aadva Regular	55066	58
External, Coral 5.1	55013	3	Hexalobe ceramic	55076	68			
Internal, Coral & Ocean 3.5, 4.1	55051	43	Conelog			Geass		
Conical, Ocean 3.5	55066	58	Screw-Line 3.3	55066	58	Way Extra	55081	73
Conical, Ocean 4.1	55004	18	Screw-Line 3.8, 4.3	55066	58			
			Screw-Line 5.0	55075	67	Global D		
B&B Implant						In-Kone 4.0-5.0	55024	15
Duravit Slim	55021	11						
Duravit 3P	55035	27						
Duravit Evolution	55035	27						
			BEGO					
Semados Mini	55012	2	Semados S / RI 3.25-5.5	55034	26			
Semados S / RI 3.25-5.5	55034	26	Semados RS / RSX 3.0	55082	74			
Semados RS / RSX 3.0	55082	74	Semados RS / RSX 3.75-5.5	55034	26			
Semados RS / RSX 3.75-5.5	55034	26	Semados SC / SCX 3.25-5.5	55034	26			
Semados SC / SCX 3.25-5.5	55034	26						

	No	Type		No	Type		No	Type
gmi			Implant Swiss			Medentis		
avantgard RP 3.75-4.25	55066	58	Tissue level	55014	4	ICX Premium 3.75 - 4.8	55062	54
avantgard WP 4.75-5.75	55075	67	Bone level Regular	55037	29	ICX Plus 3.45	55039	31
frontier RP,WP	55051	43	Bone level 3.3	55065	57	ICX Narrow	55058	50
phoenix NP 3.3	55038	30						
phoenix Std 3.3 - Wide 5.0	55011	1	Implantvision			Medical Instinct		
Monolith 3.0	55043	35	Internal Hex 3.7-5.7	55034	26	Hex 3.0	55071	63
						Hex 3.4, 4.0, 5.0	55052	44
						Cone 3.0, 3.4, 4.0	55071	63
						Cone 5.0	55072	64
GP Implants			Innobiosurg					
SSI Spiral Shape implant	55051	43	MagiCore Regular	55018	8			
gt medical			Intoss Anchor			Medigma Biomedical		
HXE	55011	1	Internal Hex 3.7-5.7	55088	80	Mars 3,5-6,0	55057	49
HXE RP	55012	2				FixTite 3,5-6	55057	49
HXE WP	55013	3				Fix-a-Dent 3,75-6	55057	49
HXI 4.1	55037	29	Intra-lock	55051	43			
002 HXI 4.7	55049	41	Intrahex 3.4 - 4.5	55051	43	Medimecca		
OCT 4.3	55014	4	Gold & Blue 3.4 - 6.0	55015	5	PT SLA Chaorum	55047	39
			InDex 3.4 - 3.75	55020	10			
			InDex 4.0 - 4.3	55012	2			
			FlatTop 3.75			Megagen		
Hahn						Anyridge	55051	43
Tapered 3.0	55058	50				Anyone Internal	55037	29
Tapered 3.5	55046	38	IRES			Mini	55064	56
Tapered 4.3	55046	38	External Hex Narrow 3.3	55012	2	MUA	55080	72
Tapered 5.0	55037	29	External Hex Regular 3.7 - 5.2	55011	1			
Tapered 7.0	55048	40	Internal Hex Narrow 3.3	55087	79			
			Internal Hex Regular 3.7 - 5.2	55049	41	Microdent		
			Internal Octagon Regular 3.7 - 5.2	55014	4	Universal	55011	1
Hiossen						Microdent 2.8 - 3.25	55068	60
HS System	55028	20	J Dental Care			Microdent 3.3 - 5.0	55011	1
ET System	55037	29	JD Evolution 3.7-6.0	55067	59	Microdent 5.5 - 8.0	55013	3
ET 3.2	55047	39	JD Evolution S 3.25	55044	36	Genius 3.0	55025	16
ET mini	55015	5	JD ICON ultra, S 2.75	55064	56	Genius 3.5 - 5.0	55051	43
			Conical Abutment	55080	72	Ektos 3.2	55023	14
						Ektos 3.7 - 8.0	55051	43
Hi-Tec			Keystone			Trylogic 3.5 - 4.2	55051	43
Tapered Self Thread NTI 3.3, SLT 3,75	55034	26	Prima / Connex Narrow	55049	41	Trylogic 5.0	55075	67
Tapered Self Thread GTI 4.2, WTI 5.0	55051	43	Prima / Connex Regular	55067	59			
Spark TPI	55034	26	Prima / Connex Wide	55041	33	MIS		
Logic Plus	55034	26	Prima / Genesis Narrow	55049	41	Conical connection Narrow (C1, V3)	55081	73
Logic LGI 3.0-3.5	55046	38	Prima / Genesis Regular	55067	59	Conical connection Standard (C1, V3)	55087	79
Logic LGI 4.3-5.0	55037	29	Prima / Genesis Wide	55041	33	Conical connection Wide (C1)	55041	33
			Restore Narrow	55079	71	Seven and M4 Narrow	55027	19
IDI Evolution			Restore Regular	55011	1	Seven and M4 Standard	55049	41
Standard implant 3.75-4.0	55004	18	Restore Wide	55013	3	Seven Wide 5.0	55051	43
B1One High	55080	72				Connect abutment NP, SP, WP	55115	81
B1One Low	55033	25	Klockner					
IDI Implants			VEGA Micro 3.0	55074	66	Neo Biotech		
ID Cam ST 4.2	55010	13	VEGA Narrow 3.5	55065	57	IS II 4.0	55020	10
ID All	55010	13	VEGA Regular 4.0, 4.5	55034	26	IS III 4.0	55020	10
			KL Narrow 3.3	55012	2			
ids			KL Regular 3.7, 4.2	55011	1	Neodent		
Reflect Aspire	55081	73	KL Wide 4.7	55013	3	Drive CM	55025	16
Reflect Certus	55019	9	SK2	55038	30	Titamax WS Cortical	55056	48
Reflect Tapered	55034	26	Essential Cone	55067	59	GM	55066	58
Reflect Recover	55020	10	S3M 3.1	55079	71			
Reflect Rapid	55020	10				Neoss		
			Konmet			SP (Ø 3.5 implants and wider)	55004	18
Implance			Classic 4.0, 4.8	55014	4	NP (Ø 3.25)	55027	19
Tissue level	55014	4	Cylindrical 3.5-4.0	55065	57	Access Abutment	55060	52
Bone level Regular	55037	29	Cylindrical 4.5-5.5	55057	49			
Bone level 3.3	55065	57	Short Cylindrical 5.0-6.0	55078	70	Nobel Biocare		
			Conical 4.4/3.4	55065	57	Brånemark System Ext hex RP	55011	1
Implant Direct			Conical 4.8/3.8, 5.5/4.5	55057	49	Brånemark System Ext hex NP	55012	2
Legacy 1/2/3/4 3.0/3.2	55012	2				Brånemark System Ext hex WP	55013	3
Legacy 1/2/3/4 3.5/3.7 - 5.7/7.0	55051	43	Lasak			Tri-lobe Connection NP	55022	12
RePlant 3.5	55022	12	BioniQ QR	55066	58	Tri-lobe Connection RP, WP	55010	13
RePlant 4.3 - 5.0	55010	13	BioniQ QN	55021	11	Tri-lobe Connection 6.0	55032	24
RePlant 6.0	55032	24				Conical Connection 3.0	55058	50
ReActive 3.5 Platform	55022	12	Leader			Conical Connection Narrow	55046	38
ReActive 4.3 - 5.0 Platform	55010	13	Internal Hex, 3.3 - 5.5 mm	55034	26	Conical Connection Regular	55037	29
Interactive 3.0 Platform	55046	38	Cono-In	55067	59	Conical Connection Wide 5.5	55048	40
Interactive 3.4 Platform	55037	29				MUA NP & RP	55080	72
Screwindirect 5.0 PF 3.0 - 5.7	55070	62	Maxillent			On1 NP, RP, WP	55042	34
			iRaise	55049	41			
Implant One						Noris Medical		
100 Series Std thread	55029	21	Medentika			Tuff 3.3 - 6.0	55034	26
200 Series Std thread	55029	21	Microcone NI	55021	11			
300 Series Std thread (blue)	55029	21	Microcone RI	55081	73	Novem		
400 Series Std thread (purple)	55035	27	Quattrocone RI	55081	73	3.4	55039	31
500 series Std thread (green)	55035	27	Quattrocone 30 AI	55081	73	3.8	55016	6
300 Series Wide thread (blue)	55029	21	Procone 3.3	55119	85	4.2, 5.0	55057	49
400 Series Wide thread (purple)	55035	27	Procone 3.8	55031	23			
500 Series Wide thread (green)	55035	27	Procone 4.3	55031	23			
			Procone 5.0	55032	24			

	No	Type		No	Type		No	Type
NucleOss			Quest			Surgikor		
Tpure Bone Level	55004	18	Hylock NP	55030	22	Fixation 3.0	55058	50
T6 Bonel Level	55074	66	Hylock RP	55066	58	Fixation 3.5	55066	58
						Fixation 4.3	55037	29
OCO Biomedical			Resista			Solution 3.5, 4.0	55066	58
Engage 3.25	55027	19	WIS 3.0-5.0	55049	41	Solution 4.5, 6.0	55037	29
Engage 4.0-6.0	55051	43				Immediate, all	55049	41
ERI 3.25	55069	61	Rex Implants			Versatile, all	55049	41
ERI 4.0-5.0	55014	4	External hex TL 1.8	55060	52			
TSI 3.25	55068	60	Internal hex TL 2.9	55034	26	Sweden & Martina		
TSI 4.0-5.0	55011	1				Syra, Syra SP 3.9-5.0	55011	1
			Ritter			Outlink ² 3.75, 4.10, 5.0	55011	1
OsseoFuse			Spiral SB/LA Narrow PF 2.9	55081	73	Shelta 3.8-5.0	55025	16
Hexaplus, Sinusplus	55081	73	Spiral SB/LA Standard PF 3.75	55049	41	Premium Straight & SP 3.3-3.8	55025	16
						Kohno Straight & SP 3.3-3.8	55025	16
Osstem			RS Implant			Premium Straight & SP 4.25-6.0	55004	18
US Regular	55011	1	3.5 - 6.0	55020	10	Kohno Straight & SP 4.25-6.0	55004	18
US Mini	55012	2				Prama 3.8 - 5.0	55079	71
US Wide	55013	3	Schutz Dental			DAT Connection	55035	27
US Wide PS	55011	1	Impla cylindrical	55052	44			
SS Regular, Wide	55028	20	IMPLA micro Retention Cone 4.2	55052	44	Tag Dental		
TS Regular, Reg Short, Diam 6 & 7	55037	29				Axis	55051	43
TS Mini 3.5	55015	5	SGS Implant			massif	55051	43
TS Mini 3.0	55047	39	P1, P7 Platform	55049	41			
			P1D, P7D	55049	41	Tatum		
Osteo Plus			P7N	55059	51	Integrity 3.7 - 5.0	55015	5
Iridium PF 3.8 & 4.2	55079	71				Integrity 6.0 - 8.0	55040	32
She 3.45, 3.75	55011	1	SIC invent			T implant	55063	55
Shi 3.75, 4.2	55049	41	SICace 3.4 - 5.0	55052	44			
Shi 5.0	55051	43	SICmax, SICtapered 3.7 - 5.2	55052	44	TBR		
			SICvantage [®] max, SICvantage tapered 3.0	55026	17	Infinity PF 3.5, PF 4.0	55004	18
OT-Medical			SICvantage [®] max, SICvantage tapered 3.7	55116	82	M Implant	55044	36
OT-F2 3.4, 3.8	55053	45	SICvantage [®] max, SICvantage tapered 4.2-5.2	55054	46			
OT-F2 4.1, 5.0	55073	65	MUA	55120	86	Thommen		
						SPI 3.5	55021	11
Overmed			Simply Implants			SPI 4.0	55065	57
Easyline Easydip	55020	10	Internal Hex	55051	43	SPI 4.5-6.0	55066	58
Oxy Implant			SIN			ticare		
Internal	55020	10	Strong SW CM	55023	14	Osseus Mini, Standard	55011	1
Kone Conical Narrow	55021	11	Strong SW Int Hex	55035	27	Osseus Maxi	55013	3
Kone Conical Regular	55051	43	Unitite CM 2.9	55084	76	Osseus Maxi ps	55011	1
PSK Line Regular	55051	43	Unitite CM 3.5	55077	69	Inhex Mini	55021	11
External Mini	55079	71	Unitite CM Compact 4.0	55085	77	Inhex Standard	55046	38
External Regular	55011	1	Tryon 3.25 - 5.0	55011	1	Inhex Maxi	55017	7
External Maxi	55011	1				Osseus Quattro Standard	55011	1
			Southern Implants			Inhex Quattro Standard	55046	38
Paltop			Tri-Nex 4.3-5.0	55010	13			
Standard 3.75-5.0	55034	26	Tri-Nex 3.5	55022	12	TRATE		
Narrow	55027	19	External Hex	55011	1	ROOTFORM	55065	57
PCA Conical	55066	58	MAX Ext Hex	55011	1			
			Co-Axis Ext Hex	55011	1	TRI Dental implants		
Phibo			Piccolo Ext Hex	55050	42	Tissue level Octa	55014	4
Aurea EVO NP	55027	19	IT-Connection	55014	4	Tri-Vent, all	55034	26
Aurea EVO RP, WP	55051	43	Deep Conical 3.0	55064	56	Narrow	55025	16
TSA S3	55074	66	Deep Conical 3.5-4.0	55046	38			
TSA S4, S5	55034	26	Deep Conical Tapered 5.0	55017	7	Trinon		
TSH S3, S4, S5	55011	1				Q2-Implant	55011	1
			Spiral Tech			QK-Implant	55014	4
P-I Brånemark			ESi, Hex connection	55051	43	QZA 3.5 Platform	55049	41
HEX 3.3	55115	81	Ultimate, Hex connection	55051	43	QZA 4.5, 5.7 Platform	55051	43
HEX 3.75 -5.0	55117	83	ESi Conical 3.5	55044	36			
AMP 3.3	55039	31	ESi Conical 4.3, 5.0	55037	29	URIS		
AMP 3.75 - 4.8	55118	84				Omni Narrow 3.0-3.5	55115	81
MT 3.3 - 4.8	55046	38	Straumann			Omni Regular 4.0-6.5	55020	10
			Tissue level, Regular & Wide Neck	55014	4			
Prodent Italia			Tissue level, Narrow Neck Crossfit (NNC)	55076	68	Vulkan		
Twinner & Prime, Conomet, Collar: Orange	55027	19	Bone level / Bone Level Tapered Narrow	55061	53	Internal Hex	55034	26
Twinner & Prime, Conomet, Collar: Fuchsia	55051	43	Bone level / Bone Level Tapered Regular	55062	54	External Hex NP	55012	2
Twinner & Prime, Conomet, Collar: Green	55067	59	Bone level / Bone Level Tapered SC 2.9	55064	56	External Hex RP	55011	1
Twinner & Prime, Conomet, Collar: Yellow	55020	10	Bone Level X (BLX) all platform	55027	19	Conical con MiP	55058	50
Prime, Conomet, Collar: Blue	55020	10	Screw retained Abutment NP, RP	55080	72	Conical con NP	55046	38
			Versatile, all	55049	41			

	No	Type
Winsix		
Internal Hex	55020	10
Xgate		
Conical 3.3 - 5.0	55034	26
Xive		
S-Plus 3.0	55030	22
S-Plus 3.4	55061	53
S-Plus 3.8-5.5	55062	54
TG Plus 3.4-4.5	55018	8
Zeramex		
XT Implant 4.2, 5.5	55086	78
Zimmer		
Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41
Screw-Vent, Tapered Screw-Vent 4.7-5.7	55051	43
Tapered Swiss Platf 3.8	55070	62
Tapered Swiss Platf 4.8	55014	4
Spline Twist Platf 4 & 5	55038	30
Eztetic NP 3.1	55065	57
Zuga		
All	55027	19