




SCALING

Model	Max
 G1	0- G 7
 G2	0- G 7
 G3	0- G 7
 G4	0- G 5
 G5	0- G 6
 G6	0- G 8
 G8	0- G 7
 G9	0- G 5
 G11	0- G 6
 G12	0- G 7
 G13	0- G 6
 G16	0- G 8

ENDODONTICS

Model	Max
 E4	0- E 6
 E4D	0- E 7
 E5	0- E 10
 E6	0- E 10
 E7	0- E 6
 E7D	0- E 6
 E8	0- E 2
 E8D	0- E 2
 E15D	0- E 6
 E11	0- E 4
 E12	0- E 3

PERIO

Model	Max
 P1	0- P 6
 P1D	0- P 10
 P2D	0- P 8
 P3D	0- P 8
 P10	0- P 10
 P11R/L	0- P 10
 P20	0- P 8
 P21R/L	0- P 8
 P25R/L	0- P 10
 P26R/L	0- P 6
 P40/P41	0- P 3
 V-Tip	0- P 6
V-P10 V-P11R/L V-P12 V-P26R/L	